

Technical Data Sheet

Article-No.: **02201**
Description: **MED-COMFORT**
facemask
3-ply, with elastic-band



SIZES	
length	width
175 mm	95 mm

PRODUCT DESCRIPTION							
material	☒ non-woven PP	☐ paper					
colour	☒ blue [B]	☐ berry [Berry]	☒ bordeaux [BR]	☒ cedro [C]	☒ green [G]	☐ lagoon [Lagoon]	☒ purple [L]
	☒ lime [Lime]	☒ Mix (blue,yellow,	☒ orange [O]	☒ pink [P]	☒ white [W]	☒ yellow [Y]	☐ platinum / grey [PT]
	☐ blue patterned [DB]	☐ wine-red patterned [WR]	☐ black [SCH]	☐ black-red-gold	☐ clean ocean / ocean green [CO]		
characteristics	☒ type IIR	☒ 3 layers	☒ latex-free	☐ with valve	☐ anti-fog	☒ with ear loops	☐ with activated carbon filter
characteristics	bacterial filtration capacity $\geq 98\%$, splash resistance ≥ 16 kPa						

REGULATORY AFFAIRS							
PPE-Regulation (EU) 2016/425	☐ category I	☐ category II	☐ category III	☒ no PPE-article			CE
MD-Regulation (EU) 2017/745	☒ class I	☐ class II	☐ class III	☐ sterile	☐ measuring fundtion	☐ no medical device	

STANDARDS							
tested in accordance to	☒ EN 14683:2019	☐ EN 149:2001					

Technical Data Sheet

Article-No.: **02201**
Description: **MED-COMFORT**
facemask
3-ply, with elastic-band



LOGISTIC DATA SUBPACKING		
general information		
material	carton	
pieces per subpacking	50	
GTIN subpacking blue [B]	4044941030753	
GTIN subpacking berry [B]	-	
GTIN subpacking bordeaux [BR]	4044941711874	
GTIN subpacking cedro [C]	4044941706948	
GTIN subpacking green [G]	4044941030760	
GTIN subpacking lagoon [Lagoon]	-	
GTIN subpacking purple [L]	4044941031019	
GTIN subpacking lime [L]	4044941029986	
GTIN subpacking mix [Mix]	4044941030982	
GTIN subpacking orange [O]	4044941030777	
GTIN subpacking pink [P]	4044941030784	
GTIN subpacking white [W]	4044941030746	
GTIN subpacking yellow [Y]	4044941030791	
GTIN subpacking platinum [PT]	-	
GTIN subpacking blue patterned [DB]	-	
GTIN subpacking wine-red patterned [WR]	-	
GTIN subpacking black [SCH]	-	
GTIN subpacking black-red-gold	-	
GTIN subpacking clean ocean [CO]	-	
measures & sizes		
length	180 mm	
width	100 mm	
height	85 mm	
weights		
size	net weight	brut weight
uni	165 g	200 g

LOGISTIC DATA PALETTE		
general information		
kind of palette	euro-palette	
measures & sizes		
cartons per layer	8	
layer per palette	8	
palette height	167	
weights		
size	net weight	brut weight
uni	160 g	185 g

LOGISTIC DATA OUTER PACKING		
general information		
material	carton	
subpackings per outer packing	10	
GTIN outer packing blue [B]	4044941030814	
GTIN outer packing berry [B]	-	
GTIN outer packing bordeaux [BR]	4044941711881	
GTIN outer packing cedro [C]	4044941706955	
GTIN outer packing green [G]	4044941030821	
GTIN outer packing lagoon [Lagoon]	-	
GTIN outer packing purple [L]	4044941031026	
GTIN outer packing lime [L]	4044941029993	
GTIN outer packing mix [Mix]	4044941030999	
GTIN outer packing orange [O]	4044941030838	
GTIN outer packing pink [P]	4044941030845	
GTIN outer packing white [W]	4044941030807	
GTIN outer packing yellow [Y]	4044941030852	
GTIN outer packing platinum [PT]	-	
GTIN outer packing blue patterned [DB]	-	
GTIN outer packing wine-red patterned [WR]	-	
GTIN outer packing black [SCH]	-	
GTIN outer packing black-red-gold	-	
GTIN outer packing clean ocean [CO]	-	
measures & sizes		
length	520 mm	
width	190 mm	
height	190 mm	
weights		
size	net weight	brut weight
uni	2.000 g	2.500 g

Technical Data Sheet

Article-No.: **02201**
Description: **MED-COMFORT**
facemask
3-ply, with elastic-band



WARNINGS AND SAFETY INFORMATION

storage / expiry date	Store mask in original packaging in a cool and dry place without additional weight, protect from direct sunlight. Do not store near ozone sources (laser printers, copiers). The actual expiry time in use cannot be specified in general terms, as it depends on the general conditions of use. An individual risk assessment must be carried out in each case. The expiry date - valid for proper storage - is stated on the packaging.
use and control	Always use face masks only for their intended use. A check/risk assessment must be carried out to ensure that the face mask is suitable for the intended use. Change the mask when it is soaked through.
disposal	Used face masks must be disposed of in accordance with the disposal regulations of the local waste disposal company. Unused face masks can be disposed of with household waste.
warnings	Change the mask when it is soaked through. Article is highly flammable - keep away from heat and/or fire.
donning and doffing instructions	PUTTING ON: Place the mask over the mouth and nose on the face and secure it in place with straps or ear loops. Adjust the flexible nose clip to the shape of the face so that as little air as possible can escape from the upper end. PULL OUT: Loosen the straps to remove the face mask from the face. Do not touch the fleece material. Disinfect your hands after disposing of the face mask.